



Cell-based Models for the Discovery and Development of Cancer Therapeutics



Utsav Sharma, PhD
Product Manager – Oncology

Hyeyoun Chang, PhD
Scientist, Cell Biology R&D

James Clinton, PhD
Lead Scientist, Cell Biology R&D

Credible Leads to Incredible™

“This Exhibitor Spotlight Theater is a promotional activity and is not approved for continuing education credit. The content of this Exhibitor Spotlight Theater and opinions expressed by presenters are those of the sponsor or presenter and are not of the American Association for Cancer Research (AACR).”




ATCC – Life Science Innovations That Touch People

A trusted partner to the global scientific community since 1925



- Founded in 1925, ATCC is a non-profit organization with HQ in Manassas, VA, and R&D & Services center in Gaithersburg, MD
- World's premiere biological materials resource and standards development organization
- World's largest, most diverse biological materials and information resource for cell biology – the “gold standard”
- Innovative R&D company featuring advanced cell models and immuno-oncology tools
- Leading global supplier of authenticated cell line, viral and microorganism standards
- Supports the global scientific community, with sales and distribution over 140 countries, 18 International distributors

Agenda



- Evolution of Cancer Therapeutics
- ATCC Oncology Portfolio
- Immuno-Oncology Reporter Models
- Patient-derived Models from the HCM1
- Q&A

ATCC credible cell products are critical to the discovery, invention, and development of incredible scientific advances

Evolution of the ATCC Portfolio




Utsav Sharma, PhD





Product Manager – Oncology, ATCC

Dr. Utsav Sharma is the Product Manager for the Oncology portfolio at ATCC. Utsav obtained his PhD in Cancer Biology at the University of Miami School of Medicine and completed his postdoctoral training at Georgetown University. His prior research work was focused on metastasis, liquid biopsy, and clinical cancer research. Prior to ATCC, Utsav worked as a Senior Scientist at Autolus Therapeutics PLC on their CD19 CAR-T platform Obe-cel, providing technical and scientific oversight to their global phase-II clinical trial on adults with Acute Lymphoblastic Leukemia (ALL) disease. In his current role, Utsav oversees all aspects of the product strategy for the oncology segment at ATCC.

Evolution of Oncology Therapies and Cell-Based Models



Oncology Portfolio



Classical Cell Culture

- ✓ Human Cell Lines
- ✓ Animal Cell Lines
 - ✓ Certified Reference Material
- ✓ Cell Culture Media
 - ✓ Cell Culture Reagents

Biomarker Discovery

- ✓ Matched Tumor - Normal Cells
 - ✓ Cancer Cell Panels
- ✓ Quantitative Cell Line DNA

Tumor Biology

- ✓ Cell Lines by Gene Mutation
 - ✓ EMT/MET Reporter Cells
 - ✓ Fluorescent Reporter Labeled Cells
- ✓ Luciferase Labeled Cells
 - ✓ Exosomes

Drug Screening

- ✓ Isogenic Cell Lines
 - ✓ EMT/MET Reporter Cells
- ✓ Primary Cells
 - ✓ hTERT-immortalized Primary Cells
- ✓ iPSC-derived Cells

Immuno-Oncology

- ✓ Primary Immune Cells and Cell Lines
 - ✓ THP-1 Reporter Cells
- ✓ Hybridoma Cells
 - ✓ iPSC-Derived Immune Cells
 - ✓ CAR-T Target Reporters
 - ✓ Checkpoint Reporter Cells
- ✓ Assay Ready Immune Cells

Patient-Derived Models


- ✓ HCMI Organoids
 - ✓ HCMI Adherent and Suspension Cell Models
 - ✓ Conditionally Reprogrammed (CRC) Cells
- ✓ Organoid Growth Kits

ATCC EXPERT
 Hyeyoun Chang, PhD
 Scientist
 R&D

ATCC EXPERT
 James Clinton, PhD
 Lead Scientist
 R&D




ATCC Immuno-oncology R&D Program



Hyeyoun Chang, PhD
Scientist, ATCC


Hyeyoun Chang, PhD, is a Scientist in the Immuno-oncology group of the R&D department at ATCC. She has extensive experience in the fields of biomedical engineering and cancer biology that focuses on drug delivery, intracellular signaling, and gene therapy. Prior to joining ATCC, Dr. Chang received her PhD in biomedical engineering from Korea University of Science and Technology and completed her postdoctoral training at Dana-Farber Cancer Institute/ Harvard Medical School.




ATCC Immuno-oncology Models

ATCC Immuno-oncology Tools

- **Develop & produce:**
 - Novel cell lines
 - Cell assays
- **Generate Application Data**
 - Conduct scientific research
 - Publish our work
- **Collaborate** with scientists from:
 - Academia
 - Small biotech/startup
 - Medium-to-large pharma




ATCC Immuno-oncology Tools



Primary and iPSC-Derived Primary Immune Cells




- Cells from most hematopoietic lineages available
- Many can be further differentiated
- Donor information available: Gender, ethnicity, age, cause of death
- iPSC-derived primary cells provide almost unlimited source of cells from same lot

Cell Type	ATCC® No.	Positive Biomarkers
CD14+ Monocytes	PCS-800-010™	CD14, CD45
Peripheral Blood Mononuclear Cells	PCS-800-011™	CD45
Bone Marrow CD34+ Cells	PCS-800-012™	CD34, CD45
Bone Marrow Mononuclear Cells	PCS-800-013™	CD45
Cord Blood CD34+ Cells	PCS-800-014™	CD34, CD45
CD4+ Helper T Cells	PCS-800-016™	CD3, CD4, CD45
CD8+ Cytotoxic T Cells	PCS-800-017™	CD3, CD8, CD45
CD19+ B Cells	PCS-800-018™	CD20, CD45
CD56+ Natural Killer Cells	PCS-800-019™	CD45, CD56
iPSC-derived MSCs	ACS-7010™	CD29, CD44, CD73, etc.
iPSC-derived CD34+ Cells	ACS-7020™	CD34, CD45
iPSC-derived Monocytes	ACS-7030™	CD14




iPSC-derived CD34+ Differentiation Potential

Macrophage differentiation capability





ATCC Immuno-oncology Tools




THP-1 Luciferase Reporter Cells

- Originate from a human monocyte-like line
- Naturally expresses pattern-recognition and cytokine receptors
- Luciferase reporter provides a straightforward, means to measure biological processes through bioluminescence measurements
- Allows broad application in:
 - Drug development,
 - Cell signaling pathway research
 - Safety evaluation for new drugs

Designation	ATCC® No.	Signaling Pathway	Function
THP-1 NF-κB-Luc2	TIB-202-NF-κB-Luc2™	NF-κB	Pivotal mediator of inflammatory response
THP-1 GAS-Luc2	TIB-202-GAS-Luc2™	JAK-STAT (Type II)	Initiates immune cell activation and recruitment
THP-1 AP-1-Luc2	TIB-202-AP-1-Luc2™	MAPK/ERK	Regulates innate and adaptive immune response
THP-1 CRE-Luc2	TIB-202-CRE-Luc2™	cAMP/PKA	Inflammatory mediator/phagocytosis modulator
THP-1 ISRE-Luc2	TIB-202-ISRE-Luc2™	JAK-STAT (Type I)	Initiates immune cell activation and recruitment
THP-1 NFAT-Luc2	TIB-202-NFAT-Luc2™	Calcineurin-NFAT	Mediates adaptive T cell activation






ATCC Immuno-oncology Tools



CAR-T Target Reporters


Chimeric Antigen Receptor T (CAR-T) Cells




Luciferase Cell Line	ATCC® No.	Tissue/Disease	Target
WIL2-S-Luc2	CRL-8885-Luc2™	B Lymphoblastoid cell	CD19
Raji-Luc2	CCL-86-Luc2™	Burkitt's lymphoma	CD19
Daudi-Luc2	CCL-213-Luc2™	Burkitt's lymphoma	CD20
Farage-Luc2	CRL-2630-Luc2™	Lymphoma	CD20
BT-474-Luc2	HTB-20-Luc2™	Breast ductal carcinoma	HER2

CAR-T Target Reporters


WIL2-S-Luc2 + CAR-T cells 24 hrs






Raji-Luc2 + CAR-T cells 24 hours



Daudi-Luc2 + CAR-T cells 24 hours




BT-474-Luc2 + CAR-T cells 24 hours





©2023 American Type Culture Collection. The ATCC trademark and trade name, and any other trademarks listed in this publication are trademarks owned by the American Type Culture Collection unless indicated otherwise.

ATCC Immuno-oncology Tools



Checkpoint Luciferase Reporter Cells




Created with BioRender.com

Qin S, Xu L, Yi M, Yu S, Wu K, Luo S. Novel immune checkpoint targets: moving beyond PD-1 and CTLA-4. *Mol Cancer*. 2019 Nov 6;18(1):155. doi: 10.1186/s12943-019-1091-2. PMID: 31690319


Designation	ATCC® No.	Biomarker	Status
HCC827-GAS-Luc2	CRL-2868-GAS-LUC2™	PD-L1	Available
MG-63-GAS-Luc2	CRL-1427-GAS-LUC2™	CD-155	Available
NCI-H1650-GAS-Luc2	CRL-5883-GAS-LUC2™	B7-H3	Available
SUP-T1 [VB]-NFAT-Luc2	CRL-1942-NFAT-LUC2™	PD-1	Coming soon
Loucy-NFAT-Luc2	CRL-2629-NFAT-LUC2™	CD96/TACTILE	Coming Soon

Checkpoint Luciferase Reporter Cells


Luciferase signal
HCC827-GAS-Luc2




Luciferase signal
MG-63 GAS-Luc2




Luciferase signal
NCI-H1650 GAS-Luc2



Luciferase signal 24h after co-culture
HCC827 GAS-Luc2 : CD8+ T = 1:1




Luciferase signal
HCC827-GAS-Luc2




T cell-conditioned media

Luciferase signal
MG-63 GAS-Luc2




T cell-conditioned media

Luciferase signal
NCI-H1650 GAS-Luc2




T cell-conditioned media

Luciferase signal 2h after co-culture
HCC827 GAS-Luc2 : CD8+ T = 1:10



ATCC Immuno-oncology Tools




ATCC Luciferase Labeled Cell Lines

- We offer a growing portfolio of commonly used tumorigenic human or mouse tumor cell lines with fluorescent or luciferase labels
- Developed by single cell cloning
- Extensively validated for
 - Cell growth
 - Stable and robust expression of reporter
- Reporter cell lines provide new versatile tools for in vitro luminescent assays and in vivo live animal bioluminescent imaging (IVIS)


Parental cell line	ATCC® No.	Species	Tissue/disease	Photons/cell/sec*	In vivo model
B16-F1	CRL-6323™	Mouse	Melanoma	336,669	Syngeneic
B16-F10	CRL-6475™	Mouse	Melanoma	207,727	Syngeneic
4T1	CRL-2539™	Mouse	Breast cancer	2,564,375	Syngeneic
EL4	TIB-39™	Mouse	Lymphoma	277,331	Syngeneic
LL/2 (LLC1)	CRL-1642™	Mouse	Lung cancer	29,104	Syngeneic
A549	CCL-185™	Human	Lung cancer	22,325,000	Xenograft
EML4-ALK A549	CCL-185IG™	Human	Lung cancer	19,356,250	Xenograft
A375	CRL-1619™	Human	Melanoma	491,308	Xenograft
KRAS G13D A375	CRL-1619IG-1™	Human	Melanoma	347,619	Xenograft
NRAS Q61K A375	CRL-1619IG-2™	Human	Melanoma	1,192,822	Xenograft
A-431	CRL-1555™	Human	Epidermoid carcinoma	603,085	Xenograft
TF-1	CRL-2003™	Human	Leukemia	112,759	Xenograft
IDH2 mutant TF-1	CRL-2003IG™	Human	Leukemia	64,768	Xenograft
HL-60	CCL-240™	Human	Leukemia	123,330	Xenograft
U-87 MG	HTB-14™	Human	Glioblastoma	8,601,250	Xenograft
IDH1 mutant U-87MG	HTB-14IG™	Human	Glioblastoma	8,165,625	Xenograft
HT-1080	CCL-121™	Human	Fibrosarcoma	2,125,000	Xenograft
HCT 116	CCL-247™	Human	Colon cancer	2,770,000	Xenograft
PANC-1	CRL-1469™	Human	Pancreatic cancer	2,175,533	Xenograft
PC-3	CRL-1435™	Human	Prostate cancer	397,640	Xenograft
LNCaP clone FGC	CRL-1740™	Human	Prostate cancer	1,187,513	Xenograft

* The *in vitro* bioluminescence (photons/cell/sec) was quantified by Xenogen IVIS Spectrum (subject to imaging and cell culture condition).

A. B16-F10-Luc2 cell line - bioluminescence **B. 4T1-Luc2 cell line - bioluminescence**



C. U87-MG-Luc2 cell line - bioluminescence





A. B16-F10-Luc2 cells (CRL-6475-LUC2™) (1x10⁶) were subjected to subcutaneous injection into the dorsal region near the thigh of C57BL/6J mice. **B.** 4T1-Luc2 cells (CRL-2539-LUC2™) were injected into the mammary fat pad of BALB/cfC3H mice. **C.** U-87 MG-Luc2 (HTB-14-LUC2™) cells (4x10⁵) were subjected to intracranial injection into the brain of nude mice. Tumor growth was monitored using a Xenogen IVIS® Spectrum. In vivo bioluminescence imaging demonstrated the establishment of the model and the progression of tumors.

Summary of IO portfolio

- Direction of IO Portfolio
- Primary and iPSC-derived Primary Cells
- THP-1 Luciferase Reporter Cells
- CAR-T Target Reporter Cells
- Checkpoint Inhibitor Cell Lines
- In vivo luciferase models

www.atcc.org/immuno-oncology





Human Cancer Models Initiative (HCMI)

Patient-derived models from the HCMI



James Clinton, PhD
Lead Scientist, ATCC

James Clinton is a Lead Scientist and group leader in ATCC's department for Cell Biology R&D and Microphysiological Systems. His team focuses on primary cells and advanced, physiologically relevant culture systems using novel technologies. His Advanced Cell Models team works on the development and commercialization of in vitro models utilizing primary and induced pluripotent stem cells, as well as exploring emerging culture techniques such as 3-D organoids and co-culture approaches.

The Human Cancer Models Initiative (HCMI)

- The HCMI is an international consortium tasked with generating 1,000 new patient-derived, next-generation cancer models (NGCMs).
- NGCMs are well-characterized and associated with a variety of clinical and molecular data.
- NGCMs are derived from patients with diverse backgrounds, disease stages, genetic backgrounds, ages, sex, treatment histories, clinical tumor diagnoses and clinical tumor stage.
- NGCMs utilize advanced culture methods that permit a diverse array of tumor tissues to be derived, propagated, cryopreserved and recovered.






HCFI founders and members



HCFI Supporting centers

- Nationwide Children's Hospital
- University of North Carolina
- Information Management Systems
- Frederick National Laboratory for Cancer Research, Leidos Biomedical Research, Inc.

Next-generation cancer models



Derivation

- All from primary patient tissues

Molecular characterization

- 15X whole genome sequencing
- 150X whole exome sequencing
- 120 million read RNA-sequencing
- Infinium MethylationEPIC DNA Array

Clinical data

- Standardized case report forms
- Patient demographics
- Disease diagnosis, treatment and outcome information

ATCC: Supporting the HCMI since 2016



**Authentication
and quality
control testing**



**Manufacturing
and cell bank
generation**



**HCMI model
information
resource**



**Distribution of
models through the
ATCC catalog**

Authentication and quality control testing

Standard testing performed on all HCMI models

- STR (short-tandem repeat) analysis
- CO1 (Cytochrome Oxidase C1) testing
- Mycoplasma contamination detection
- Sterility testing (BACT/ALERT® 3D)
- Human virus panel (HIV, HepB, HPV, EBV, and CMV)



ISO 9001:2015 Certification for quality management system


- Demonstrates commitment to quality products, customer service, and continued improvement



ISO/IEC 17025:2017 Accreditation for testing

- Applies to all ATCC cultures, derivatives, and bioproducts tested in our laboratories

Manufacturing and cell bank generation



HCMC model information resource

HCMC Landing page atcc.org/hcmi

Human Cancer Model Initiative

277 Products
Results 1-24 of 277

Search

Filtering options

- Product category
 - Human cells 270
 - Media 7
- Product type
 - Cell model 270
 - Organoid 180
- Product application
 - 3D cell culture 245
 - Cancer research 235
 - Neuroscience 44
 - Cell culture 7
 - Cell growth and viabi... 7
- Product format
 - Frozen 277
- Biosafety level
 - BSL 1 258
 - No BSL 7
 - BSL 2 2
- Organism
 - Homo sapiens 270

HCMC model-specific product page

Product: HCM-BROD-0611-C71
PDM-401 BSL 1

Product format: Frozen
Tissue: Skin
Disease: Glioblastoma Recurrent

QUICK VIEW Compare

HCM-CSHL-0065-C20
PDM-9*

Product category: Human cells
Product type: Cell model
Organism: Homo sapiens, human
Morphology: Organoid
Tissue: Large Intestine, Rectum
Disease: Sarcomatoma, Primary metastatic
Applications: 3D cell culture, Cancer research
Product format: Frozen
Storage conditions: Vapor phase of liquid nitrogen

Purchase model

Price: \$3,634.00
Quantity: **ADD TO CART**

Model resources and protocols, lot-specific information

- Product sheet
- Certificates of analysis
- Safety data sheet
- Protocol
- Protocol
- Protocol

Detailed product information

- General
- Characteristic
- Handling information
- Quality control specifications
- History
- Legal disclaimers

Permits & Restrictions

Material Transfer Agreement Addendum for Organoid Products
For every order of this item, you must provide a signed [Material Transfer Agreement Addendum for Organoid Products](#). We cannot ship this item until we receive this addendum. Email the signed addendum to [customerservice@atcc.com](#) with a reference to your account and order numbers. Once received, your addendum will be reviewed, and this item will be released for shipment if all requirements are met. If you need assistance with your order, please contact our Customer Care team at your applicable distribution.

Import Permits for the State of Hawaii
If shipping to the U.S. state of Hawaii, you must provide written an Import permit or documentation stating that an Import permit is not required. We cannot ship this item until we receive the documentation. Contact the [Hawaii Department of Agriculture \(DOA\), Plant Industry Division, Plant Quarantine Permit](#) to determine if an Import permit is required.

QC testing specifications

HCMC Searchable Catalog


GDC

Model images and other information

ATCC Number: HCM-CSHL0065-C20
Designation: HCM-CSHL0065-C20


HCM1 model information resource

Validated companion reagents




- Seven HCM1 model formulation-specific reagent kits
- Extracellular matrix
- Growth factor secreting cell lines

Video protocols, other educational resources



Also, application notes, webinars, poster presentations, and troubleshooting tips.

HCMI models in the ATCC catalog




Recent releases

- First head and neck cancer model
- First desmoid tumor cell line

In the production pipeline

Additional lung, breast, endometrial, ovarian, kidney and first bladder models

HCMI portfolio model diversity







TOP 10 (80% of total)

Brain	49
Pancreas	46
Colon	40
Esophagus	32
Skin	26
Rectum	16
Stomach	13
Connective tissue	9
Biliary tract	8
Lung	6


Collection includes models derived from rare adult and pediatric cancers such as rhabdomyosarcoma, leiomyosarcoma, Ewing sarcoma and Wilms tumor.

Complete list available at
atcc.org/hcmi

HCMI portfolio clinical and demographic characteristics



Age distribution




Increasing patient diversity of models is a goal for on-going model releases.


HCMI data aligns with patient genotypes

Top 10 frequently mutated genes vs. The Cancer Genome Atlas


Colon




Skin



Brain



Pancreas



Top 10 genes affected in the HCMI affected cases

APC	LRP1B	MUC16	GRIN2A	PTEN	PIK3CA	KRAS	MUC4
TP53	ZFHX4	LRP1B	MECOM	TP53	CSMD3	TP53	MUC16
KRAS	CSMD3	CSMD3	FAT4	EGFR	BCOR	SMAD4	KMT2D
MUC16	ROBO2	FAT3	DCC	CNTNAP2	NF1	CDKN2A	FAT3
FAT4	ALK	FAM135B	COL1A1	PIK3CR1	FAT1	ACVR2A	RHOH

Also, TP53, BRAF, NRAS, NF1

Canonical mutation

Matching mutation

Help us prioritize future HCMI model releases

Hundreds of models have not yet entered the ATCC HCMI manufacturing pipeline.

Browse and search unreleased HCMI models at ATCC

- Use the “Submit your Input” button on the HCMI Landing page
 - Direct link: www.atcc.org/hcmi-input
- Look for “Expansion Status” on the HCMI Searchable Catalog

Email us which HCMI models are most relevant for your research

- Contact us at: hcmi@atcc.org



Resources to learn more about ATCC and the HCMI

The **Human Cancer Models Initiative (HCMI)** is an international consortium that is dedicated to generating novel human tumor-derived culture models with associated genomic and clinical data. The **HCMI consortium** comprises funding agencies and cancer model development institutions. The consortium's funding agencies include the **National Cancer Institute (NCI)**, **Cancer Research UK (CRUK)**, **Hubrecht Organoid Technology (HUTO)**, and **Wellcome Sanger Institute (WSI)**. NCI-funded model development institutions include the **Broad Institute** and the **Cold Spring Harbor Laboratory**; CRUK and WSI co-fund organoid development in the United Kingdom; CRUK provides the patient samples, while WSI derives and sequences the organoid models. The foundation HUTO is a Netherlands-based not-for-profit organization that derives and sequences organoid models. ATCC was selected as the sole distributor for the HCMI models. The generating institutions deposit the models into ATCC, where they are authenticated, expanded, preserved, and made available for global distribution. The **HCMI model data** are available from the NCI as a resource to the research community.

Name	Primary Site	Clinical Tumor Diagnosis	Tissue Status	Age At Acquisition (Years)	Age At Diagnosis (Years)	Has Multiple Models	Expansion Status	# Mutated Genes
HCM-BROD-0227-C43	Skin	Melanoma	Metastasis	40	40	No	EXPANDED	3075
HCM-BROD-0569-C43	Skin	Melanoma	Metastasis	79	78	No	EXPANDED	2886
HCM-CSHL-0426-C18	Colon	Colorectal cancer Primary		73	72	No	EXPANDED	2701
HCM-BROD-0097-C34	Bronchus and l.	Lung cancer	Metastasis	66	65	No	EXPANDED	2338
HCM-CSHL-0459-C17	Small intestine	Rare cancers Primary		57	57	No	EXPANDED	2426
HCM-BROD-0223-C43	Skin	Melanoma	Metastasis	74	73	No	EXPANDED	2187
HCM-BROD-0106-C71	Brain	Glioblastoma	Recurrent	56	52	No	EXPANDED	2122
HCM-BROD-0334-C43	Skin	Melanoma	Metastasis	72	70	No	EXPANDED	1619
HCM-CSHL-0174-C22	Intrahepatic bile	Intrahepatic bile Primary		64	64	No	EXPANDED	1568
HCM-CSHL-0317-C18	Colon	Colorectal cancer Primary		65	64	No	EXPANDED	1502
HCM-BROD-0025-C16	Stomach	Stomach cancer Primary		74	73	No	EXPANDED	1330
HCM-BROD-0679-C43	Skin	Melanoma	Metastasis	69	68	No	EXPANDED	765

HCMI Searchable Catalog
<https://hcmi-searchable-catalog.nci.nih.gov>

NCI Genomic Data Commons
<https://portal.gdc.cancer.gov/projects/HCMI-CMDC>

SESSION NIH03 - Human Cancer Models Initiative: A Resource of Next-generation Cancer Models and Data

Today at 2:30-3:30 PM in Room W314A




Summary and Resources



- *ATCC is the comprehensive solution provider.*
- ATCC has a vast collection of **authenticated and characterized** bio-materials intrinsic to oncology research and drug discovery.
- We fully characterizes our **cell models** at the genome, proteome, and functional levels.
- We have an established R&D program focused on the development of the next generation of advanced cell models to streamline **oncology/immuno-oncology discovery** and drug screening platforms.
- We offer **exhaustive technical and application data** along with companion products with our cell models to support your experimental workflows.

<https://www.atcc.org/coming-soon-products>



New Products:



CAR-T Target Reporter-Labeled Tumor Cells

- Access CAR-T potency and efficacy
- High endogenous expression of CAR-T target antigens
- Available for CD19, CD20, and HER2

Checkpoint Luciferase Reporter Cells

- Enables screening of checkpoint inhibitor molecules
- Wide range of targets such as PD-L1/2, CD-155, B7-H3, and PD-1
- Luciferase will be expressed under the control of GAS or NFAT

Human Cancer Models Initiative (HCMI)

- 2-D and 3-D patient-derived models available
- Diverse genetic backgrounds of the same cancer types
- Culturing protocols and organoid growth kits

