

ATCC®

Credible leads to Incredible™

The Importance of Standardizing Lyophilization Techniques in Microbial Preservation

Jyoti K. Jha, PhD,* Rahul Tevatia, PhD; Lukas Underwood, PhD; Quinn Osgood, BSE; Emmanuel Boafu-Asare, MS; Nilay Chakraborty, PhD

ATCC, Manassas, VA 20110

Abstract

Storage of microorganisms in a quantitative format is essential for modern industrial quality control and other related applications. Lyophilization is one of the most widely used techniques for preserving microorganisms in a desiccated format as it allows lab workers to avoid expensive ultra-cold chain storage and related distribution challenges. However, lyophilization is a complex process that involves phase changes of the extracellular environment, and the outcome is highly influenced by the compositional, thermal, and pressure changes of the materials. In this study, we evaluate how changes in lyophilization formulation affect microbial viability and protein expression characteristics using *Escherichia coli* (ATCC® 8739™) as a model. Our results highlight the need for optimization of lyophilization formulation for desiccated storage of microorganisms. Further studies are required to identify the standardized lyophilization processes.

Introduction

Lyophilization is a well-established method for supporting the long-term storage of bacterial strains and it is by far the most frequently used method for achieving long-term shelf life.^{1,2} In this process, the addition of a selected buffer formulation to the biomass and subsequent lyophilization are vital steps for maximizing the viability of a microorganism.³

In this study, we used *E. coli* (ATCC 8739) as model organism for evaluating viability optimization and long-term storage in a novel buffer formulation. This bacterial strain is a WT version of *E. coli*, and it has wide applications as a quality control strain in testing antimicrobial handwashing formulations, media, efficacy, and bio-resistance. Further, it is cited as a reference material in multiple United States Pharmacopoeia (USP) methods, namely USP <51> and <62>.^{4,5} Here, we used global proteomic analysis to provide evidence that may co-relate for long-term viability of *E. coli* in our novel lyophilization formulation.




Figure 1: Overview of microbial lyophilization. Lyophilization, or freeze drying, is the process of removing water from a product in its frozen state under a vacuum. This process allows ice to change from a solid to a vapor without going through a liquid phase. This process comprises three major steps: freezing the sample, primary drying (ice sublimation), and secondary drying (unfrozen water desorption).

Materials and Methods

Lyophilization Formulation for *E. coli* Preservation: Initial viability screening was performed using a combination of multiple sugars, polyols, proteins. The viability of *E. coli* preserved in buffer #1 (B#1; our optimized mixture of carbohydrates, polyols, and proteins) was compared to a commonly used lyophilization buffer, buffer #2 (B#2; a mixture of carbohydrate and proteins).

Bacterial Viability Assay: Bacterial viability was determined using a plate-based dilution method. Viable bacteria were quantitated via colony formation units (CFU).

Quantitative Proteomics Analysis (TMT-Multiplex Labeling and LC-MS/MS): To identify potential protein profile change related to different buffer and storage conditions, a comparative proteomic analysis was performed. A total of 12 samples from four conditions were analyzed in triplicate using the TMT-13plex labeling-based quantitative proteomics method and liquid chromatography-tandem mass spectrometry (LC-MS/MS). In this study, 2658 proteins were quantitatively identified. About 22 of the 2658 proteins significantly changed in abundance (expression) (≥ 2 fold; $p < 0.05$; $n=3$) in at least one condition compared to another.

Results




Figure 2: Stability of *E. coli* in two different formulations. (A-C) Stability of *E. coli* in B#1 and B#2 at 4°C, 22°C, and 37°C, respectively, for one year. Experimental data was collected for 6 months, and predictive modeling was performed to extrapolate stability up to a year. (D-F) Viability loss ($\Delta \log$ CFU) for *E. coli* in B#1 and B#2 at 4°C, 22°C, and 37°C, respectively, after 6 months of storage.




Figure 3: Sample preparation for lyophilization, viability determination, and LC-MS/MS. (A) Flow chart of sample processing for LC-MS/MS and viability of *E. coli*. (B) Viability determination for lyophilized samples of *E. coli* in B#1 and B#2 and incubation at respective temperatures (4°C and 37°C) for 6 weeks. (C) Mass spectrometry-based quantitative proteomic profiling workflow (lysate preparation, trypsin digestion, TMT-13plex labeling, purification, and analysis).




Figure 4: Proteomic characterization of *E. coli* in different storage condition and formulation. (A) A heat map showing the relative abundance of 2658 proteins identified across 12 samples from 4 group conditions depicting a total of 122 changed proteins (fold change $\geq 25\%$ ($p < 0.05$, $n=3$)). (B) Volcano plot demonstrating the fold change of protein abundance. [Note that the sample comparison in (A) and (B) is among B#1@4°C and B#2@4°C, B#1@4°C and B#1@37°C, B#1@37°C and B#2@37°C, or B#2@4°C and B#2@37°C.] (C) Select changed proteins between B#1@4°C and B#2@4°C (> 2 fold). (D) Select changed proteins between B#1@37°C and B#2@37°C (> 2 fold). Orange highlighted proteins are common changes in B#2 vs B#1 at both 4°C and 37°C. Change of other proteins are related to their conditions. No protein was found changed (> 2 fold) within the same formulation (B#1@4°C vs B#1@37°C and B#2@4°C vs B#2@37°C).

Conclusions

Here, we have identified a formulation for stabilizing the lyophilized product of *E. coli* at 4°C, 22°C, (room temperature), and 37°C storage. The viability of the lyophilized *E. coli* was found to be 2 to 3 log better for up to 6 months in our optimized formulation (B#1) as compared to a conventional formulation (B#2). We exploited LC-MS/MS to understand the proteomic variances that might have caused reduced viability of lyophilized *E. coli* in B#2. Interestingly, we found that the cold shock proteins were induced three times more in formulation B#2 as compared to B#1, which indicates a reduction in global protein synthesis⁶ and may point toward lower viability after lyophilization. Also, downregulation of the 'Ada' protein⁷ and transcriptional activator 'FaeR' in B#2 will lower the capability of DNA methylation repair and increase the chance of abnormal cell production. Further, multiple stress response proteins (BhsA) were overexpressed, which indicated that the *E. coli* cells in formulation B#2 are unstable. Elevated expression of the type III restriction modification enzyme (mod) in formulation B#2 indicated DNA cleavage and hence may be another factor that contributed to lower viability.⁸ Global mass spectrometry analysis indicates several biomolecular mechanisms contributing to variations in the viability of lyophilized *E. coli* in formulations B#1 and B#2. Further investigation is required to determine the direct impact of these proteomic differences on lyophilization stability.

References

1. Cárcoba R, Rodríguez A. Influence of cryoprotectants on the viability and acidifying activity of frozen and freeze-dried cells of the novel starter strain *Lactococcus lactis* ssp. *lactis* CECT 5180. *Eur Food Res Technol* 211: 433-437, 2000.
2. Carvalho AS, et al. Relevant factors for the preparation of freeze-dried lactic acid bacteria. *International Dairy Journal* 14(10) 835-847, 2004.
3. Boehm AB, et al. Faecal indicator bacteria enumeration in beach sand: a comparison study of extraction methods in medium to coarse sands. *J Appl Microbiol* 107(5): 1740-1750, 2009.
4. USP <51> Antimicrobial Effectiveness Test
5. USP <61> Microbiological Examination of Non-sterile Products: Microbial Enumeration Tests
6. Al-Fageeh MB, Smales CM. Control and regulation of the cellular responses to cold shock: the responses in yeast and mammalian systems. *Biochem J* 397(2): 247-259, 2006.
7. He C, et al. A Methylation-Dependent Electrostatic Switch Controls DNA Repair and Transcriptional Activation by *E. coli* Ada. *Mol Cell* 20(1): 117-129, 2005.
8. Mücke M, et al. DNA cleavage by type III restriction-modification enzyme EcoP151 is independent of spacer distance between two head-to-head oriented recognition sites. *JMB* 312(3): 687-698, 2001.



Order the products used in this study